


Changing a battery in your hearing aid

This leaflet is about changing the battery in your hearing aid


How to switch your hearing aids on



The battery drawer is at the bottom of the hearing aid.



1. Open the battery draw using your finger nail.



2. Take out the battery. This is shown by number 1 in the picture.



3. Put the old battery somewhere safe to be recycled.

4. Keep the hearing aid down on a table with the battery draw open.




5. Take a new battery out of the packet.


6. Remove the sticker from the new battery.



7. Leave the new battery for **5 minutes**. This will make it last longer.



8. Put the new battery in the draw of the hearing aid with the flat side facing up. This is shown by number 2 on the picture.



9. Gently close the battery door.



The hearing aid is now turned on and ready to be used.