


Putting your earmould in

This leaflet is about putting your hearing aid on and cleaning it

How to put your earmould in




The **earmould** is the part that goes in your ear.



Pick up one of your hearing aids. The left one or the right one.




1. Hold the **earmould** with your finger and thumb.



2. Put the end with the hole into your ear.



3. Get the piece at the top in.




4. Push the **earmould** firmly with your finger. You may need to pull down your **earlobe** too. Your **earlobe** is the part at the bottom of your ear. Where an earring may sometimes be.




5. Slide the hearing aid over the top of your ear.

Cleaning your hearing aid



Wipe the **earmould** with a tissue every evening.



Use the loop on the cleaning tool to remove any wax or skin from the hole.



Use the brush on the cleaning tool to remove any wax or skin from the hearing aid.



The tubing on the hearing aid will need to be changed about every 4 months.



- You can come to the **Audiology** clinic for help to change the tube.



- Or we can arrange for a volunteer to visit you at home. A volunteer is someone who helps for free.