


Switching your hearing aid on

This leaflet is about using your hearing aid

How to switch your hearing aids on



The battery drawer is at the bottom of the hearing aid.


Check that the battery is in.

Switch the aid on by closing the battery drawer fully.



Any whistling sounds should stop when the aid is in your ear correctly

When you have 2 hearing aids



A red dot means it is for the right ear.



A blue dot means it is for the left ear.



The **earmould** will only go in the correct ear. They won't fit in the wrong ear.

The **earmould** is the part that goes in your ear.


How to take your hearing aid out



Push the hearing aid off the top of your ear.

Pull the **earmould** out of your ear.

Switch the aid off by opening the battery drawer.



The battery drawer is at the bottom of the hearing aid.

You can use your fingernail to help to open it.


Changing the battery



The battery will need to be changed every 2 weeks.



You can wait until the hearing aid beeps, then change the battery.



Or you can ask someone else to change the battery on the same day every week.