


# Washing your earmould

This leaflet is about washing your earmould

## How to prepare your earmould for a clean



The earmould is the part that goes in your ear.



For a good clean, take apart the hearing aid and **earmould**.

Pull the soft tubing off from the hearing aid elbow. The part that is the elbow is shown in the picture.

Do not pull the tubing out of the **earmould**.

## How to wash your earmould


Wash the  **earmould**  part in warm soapy water.



Rinse well in clean water.

Shake to get rid of any water from the tube.

Leave it to dry overnight in a warm, dry place. Or wipe it dry with a tissue.



Put the tubing back onto the hearing aid elbow with the curves together.